




Patient ID:


Date:


	32301	32301/HX
	32301/NOLS	32301/MR
	32301-PED	32204
	32301-PED/SS	32204-PED
	32301-PED/NOLS	
	32301-PED/SS/NOLS	
AIO base plate		


	LP-KS-B (29026)	LP-KS-W (29027)
	LP-KS-E (29028)	LP-EC (29029)
Knee & leg positioning set		




	TA-WE-0 (29007)	
Wedge 0°		


	Upper frame (29053/1)	
	Middle frame (29053/3)	
	Bottom frame (29053/2)	
Indexation frame for knee & leg positioning cushions		


	TA-LS-0 (29011)	
Lateral support 0°		


	P-SE-LD (32140/12)	
	P-SE-HD (35754/6N)	
Secondary base plate		


	TA-WE-5 (29008)	
	TA-WE-10 (29009)	
	TA-WE-15 (29010)	
Wedges 5°, 10° and 15°		

	AC-GP-L (32015/9/2) Grip Pole Long	Grip pole collar (32015/9/:
	AC-GP-S (32015/9/3) Grip Pole Short	
Grip poles		
		

	TA-LS-5 (29012 + 29013)	
	TA-LS-10 (29014 + 29015)	
	TA-LS-15 (29016 + 29017)	
Lateral supports 5°, 10° and 15°		

	AC-GP-H (32015/10)	
Horizontal grip pole		

	TA-AR-H (29018)	
	TA-AR-L (29019)	
Arm rest low and high		

	AC-PG (32015/5)	
Protruding hand grip		

All Rights reserved. No part of this publication may be produced, stored in a retrieval system or transmitted, in any way, form or by any means, electronic, mechanical, photocopying or otherwise, without prior written permission of Orfit Industries ©ORFIT Industries





AIO Solution - Breast- & Lungboard Components


Individual patient components list


Patient ID:

Date:

	35755 (Model 1)	35765 (Model 2)
	35753 (Model 3)	35753ZF (Model 4)
	35752 (Model 5)	35752ZF (Model 6)
Set of low density head supports		

	35758-MD (Model 1)	35758ZF-MD (Model 2)
	35713-MD (Model 3)	35713ZF-MD (Model 4)
	35714-MD (Model 5)	35714ZF-MD (Model 6)
Set of regular density head supports		

	Upper frame (29054/2)	Bottom frame (29054/1)
	Position Knee Support:	
	Position Feet Support:	
AIO indexable knee and foot frame		

	29112	
	Arm Support Cushion	

Positioning of elements			
AIO base plate position (position of laser line on millimeter scale):			
Height of the protruding or horizontal handgrip:			
Position of the grey knee & leg positioning set:	A-position		
	B-position		
	Position Upper Frame	1 to 28	
	Position Bottom Frame	1 to 25	
Thermoplastic masks used:	4-points breast mask		
	5-points assymetric breast mask		
	6-points thorax-abdomen mask		
Others/remarks			

All Rights reserved. No part of this publication may be produced, stored in a retrieval system or transmitted, in any way, form or by any means, electronic, mechanical, photocopying or otherwise, without prior written permission of Orfit Industries ©ORFIT Industries